

BBL**PEPTONES**

Acidicase **Biosate**
Gelysate **Lactalysate**
Milk-Protein Hydrolysate
Myosate **Phytone**
Polypeptone
Thiotone **Trypticase**

•

Reprint Sent on Request

**BALTIMORE BIOLOGICAL
LABORATORY
INCORPORATED**

*A Division of Becton, Dickinson & Co.***Baltimore 18, Md.****CUMULATIVE INDEX**

to the JOURNAL

of BACTERIOLOGY

VOLUMES 31 TO 64 • \$7.50

Increase the usefulness of your volumes of the journal. Locate the exact information you want, quickly and easily. Turn to any subject, any author from any issue with a minimum of time and effort.

THE WILLIAMS & WILKINS CO.

**MT. ROYAL AND GUILFORD AVES.
BALTIMORE 2, MARYLAND**

.....

*Two Books of Interest to Every Bacteriologist***BACTERIA**

By K. A. BISSET, Ph.D., D.Sc.

A natural history of bacteria
Helps you to recognize and understand unfamiliar species when you meet them, by means of carefully-framed descriptions of the bacteria themselves, and full pictures of entire classes

Correlates morphological, biochemical and ecological characters

Emphasizes the author's thesis of the natural relationships of bacteria: that the characters and place in nature of any genus are implicit in, and predictable from, its systematic position

*131 pp., 38 figs., \$4.00**Revised 2nd Edition***The CYTOLOGY and LIFE HISTORY of BACTERIA**

By K. A. BISSET, Ph.D., D.Sc.

Presents bacteria as living cells with the same structure and function as other living cells

Correlates the available information on the various types of bacteria

"Very informative and intensely interesting."—*N. Y. State J. Med.*


*176 pp., 68 figs., \$6.00***THE WILLIAMS & WILKINS COMPANY**

Mt. Royal and Guilford Aves.

Baltimore 2, Maryland

(In writing to advertisers, please mention the journal it helps.)

International STANDARDS



SIZE 1, MODEL SBV CENTRIFUGE

VERSATILITY — accommodates

- 11 Horizontal swinging heads
- 15 Angle heads
- 2 Multispeed attachment heads
- 4 Basket style heads

SPEEDS and FORCES —

- Adjustable speeds up to 5600 r.p.m.
- R.C.F. — up to 5100 G.
- Multispeed attachment provides up to 20,000 r.p.m. and 27,500 G.

CAPACITY —

- Standard glass and plastic containers from 10 ml. to 250 ml.
- Total maximum capacity — 1000 ml.

SIZE 2, MODEL V CENTRIFUGE

VERSATILITY — accommodates

- 14 Horizontal swinging heads
- 16 Angle heads
- 2 Multispeed attachment heads
- 4 Basket style heads

SPEEDS and FORCES —

- Adjustable speeds up to 5100 r.p.m.
- R.C.F. — up to 4060 G.
- Multispeed attachment provides up to 19,000 r.p.m. and 25,000 G.

CAPACITY —

- Standard glassware from 10 ml. to 650 ml.
- Plastic containers up to 1000 ml.
- Total maximum capacity — 4 liters.

Throughout the world there are standards of comparison for many products. For Laboratory Centrifuges that standard of comparison is, invariably, International's Size 1 or Size 2. These Centrifuges are truly *International Standards*.

International Equipment Company offers you the experience of more than fifty years of specialization in the design and manufacture of Centrifuges. Motors and accessories, *specifically designed for these machines*, are completely manufactured in the International Plant, eliminating the chance of obsolescence of the whole machine for lack of one of its parts. Accessories are thoroughly tested and rated for the model in which they are to be used. To assure "follow-up" safety for the user, Speed and Force Tables are provided, giving performance expectations and safety limits for each combination of accessories.

With such built-in reliability, it is small wonder that bid specifications, government or private, the world over, call for these International Centrifuges.

International  Equipment Co.

1284 SOLDIERS FIELD ROAD • BOSTON 35, MASSACHUSETTS

In writing to advertisers, please mention the journal it helps.



PPLO

ISOLATION and CULTIVATION

THE AVAILABILITY of media for the isolation and cultivation of *Pleuropneumonia-Like Organisms* (PPLO) has increased the interest in these infectious agents. PPLO are readily grown from clinic specimens or from pure culture on Bacto-PPLO Agar and Bacto-PPLO Broths enriched with Bacto-PPLO Serum Fraction or Bacto-Ascitic Fluid.

BACTO-PPLO AGAR is an especially prepared solid medium containing the essential factors from meat infusion required by the PPLO. This medium, enriched with Bacto-PPLO Serum Fraction, will support abundant growth of all PPLO strains.

BACTO-PPLO BROTH w/o CV enriched with Bacto-Ascitic Fluid or Bacto-PPLO Serum Fraction is a liquid medium for enrichment and cultivation of PPLO organisms. Selective media for PPLO may be prepared from Bacto-PPLO Broth w/o CV by the addition of thallium acetate, potassium tellurite, penicillin or other bacteriostatic agents as desired.

BACTO-PPLO BROTH with CV enriched with Bacto-Ascitic Fluid is a selective medium for the isolation and cultivation of PPLO from specimens contaminated with commensal organisms. Contains Crystal Violet for inhibition of Gram-positives. The selectivity of the medium may be further increased to inhibit Gram-negatives by the addition of Bacto-Chapman Tellurite Solution. This medium will support growth of all PPLO strains and effectively suppress the growth of most contaminating bacterial flora.

BACTO-PPLO SERUM FRACTION contains the growth accessory factors required by all strains of the PPLO and is used as an enrichment in Bacto-PPLO Agar and Bacto-PPLO Broth w/o CV.

BACTO-ASCITIC FLUID is a sterile, bile-free ascitic fluid specially prepared for enriching culture media used in the cultivation of the PPLO group and other fastidious microorganisms.

BACTO-CHAPMAN TELLURITE SOLUTION to be added to Bacto-PPLO Broth w/CV or Bacto-PPLO Broth w/o CV for the inhibition of Gram-negative microorganisms.

*THE DIFCO MANUAL, NINTH EDITION,
including descriptions of these media and their use,
is available on request.*

DIFCO LABORATORIES

DETROIT 1, MICHIGAN